

Chambre d'inhalation, de faible volume (60 mL), pour circuit de ventilation mécanique.
Low volume (60 mL) inhalation chamber for mechanical ventilation circuit.
Inhalationskammer, kleines Volumen (60 ml), für den maschinellen Beatmungskreislauf.
Camera di inalazione, a basso volume (60 mL), per circuito di ventilazione meccanica.

PI-NT-MIN-0002 / 20240709


1




2

3

4

Fabriqué en France par / Manufactured in France by / Hergestellt in Frankreich von / Prodotto in Francia da :
 OptimHal-ProtexSom SAS
 24 rue du Train Renard - ZA d'Armanville BP46
 50700 Valognes - France
 +33 (0) 233 887 094
 info@optimhal.com
 www.optimhal.com

10°C  30°C

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INTENDED USE

MinimHal® is an inhalation chamber intended for adult or pediatric patients with invasive or non-invasive mechanical ventilation who have been prescribed aerosolized treatment from a pressurized metered-dose inhaler and/or a vibrating mesh nebulizer used simultaneously or not.

NOTES

- MinimHal® should be used by clinicians qualified in the care of mechanically ventilated patients.
- Please read the instructions for use of MinimHal®, the metered dose inhaler and the nebulizer and keep them available at all times of use.

COMPOSITION (FIG.1)

a. Male connector 15 mm
 b. Female connector 22 mm
 c. One opening for MDI
 d. One opening for vibrating mesh nebulizer
 e. Two sealed caps
 f. One nozzle
 g. A 15 mm female to 22 mm male converter that allows the use of MinimHal® on a pediatric circuit or on an adult circuit.
 h. A RespiMat® inhaler adapter can also be supplied.

REQUIREMENTS BEFORE USE

When putting the device into operation, it is imperative to ensure that:

- The device does not contain any foreign body, to avoid any risk of damage to the patient's airways or obstruction of the medical device.
- Outside of drug administration periods, ensure that the caps are properly closed to guarantee the airtightness of the device and prevent any risk of leakage.

INSTRUCTIONS FOR USE

Checking the operation of the nozzle
 Remove the MDI canister from its actuator. Open the orifice cap (c) and insert the canister (Fig.2).
 Operate the MDI and check that the nozzle correctly ensures its function (visually test with a puff of drugs) (Fig.3).

INSERTION IN THE CIRCUIT (FIG.4)

For insertion of MinimHal® on the inspiratory limb of a pediatric dual limb circuit:

The 15mm male connector must be connected to the inspiratory limb. The 22mm female connector must be connected to the converter connected to the Y-piece.

For insertion of MinimHal® on the inspiratory limb of an adult dual limb circuit:

The 15mm male connector must be connected to the converter itself connected to the inspiratory limb. The female connector must be connected to the Y-piece.

For insertion of MinimHal® between the Y-piece and the endotracheal tube of an adult or pediatric dual limb circuit:

In case of an invasive ventilation, the male connector must be connected to the Y piece and the female connector must be connected to the converter. The converter must be connected to the patient's interface. In case of a non-invasive ventilation, it is possible to connect the female connector directly to the patient's interface.

**When using MinimHal® after the Y-piece, make sure that the dead volume ratio (to which the chamber's volume is added - 60 mL) is compatible with the tidal volume delivered, and if necessary, disconnect MinimHal® from the circuit.*

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ENGLISH

DRUG ADMINISTRATION BY PRESSURIZED METEDED DOSE INHALER

Step 1: Shake the MDI canister before use in accordance with the instructions provided by the manufacturer. Remove the cap from the opening (c) and insert the canister (Fig.5).
Step 2: Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the MDI (Fig.6).
Step 3: Remove the canister from the MDI port and replace the cap on the orifice, shake the canister and repeat steps 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the aerosol dispenser for the waiting time between each dose. At the end of the drug administration, remove the canister from the opening MDI and replace the cap on the orifice (c) (Fig.7).

DRUG ADMINISTRATION BY VIBRATING MESH NEBULIZER

Remove the cap from the orifice (d) and insert the vibrating mesh nebulizer (Fig.8). To avoid an opening of the circuit, this step could be done before the insertion of MinimHal® into the circuit. Pour the drug solution into the nebulizer reservoir. For proper operation of the nebulizer, follow the instructions provided by the manufacturer. Aerogen® vibrating mesh nebulizers may remain on the circuit as recommended by the manufacturer.

DRUG ADMINISTRATION WITH THE RESPIMAT® INHALER

Step 1: Remove the cap from the orifice (d) and insert the adapter Respimat®. Insert the Respimat® inhaler into the adapter (Fig.9).
Step 2: To operate the inhaler, follow the instructions provided by the manufacturer. Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the inhaler.
Step 3: Repeat step 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the inhaler dispenser for the waiting time between each dose. At the end of the drug administration, remove the inhaler and the adapter from the opening (d) and replace the cap on the orifice (d) (Fig.10).

CLINICAL BENEFITS – DEVICE PERFORMANCE

The chamber volume allows the propellant gas of the larger particles size to evaporate, producing smaller aerosol particles of the drug. Smaller particles improve the efficiency of aerosol delivery to the patient through the ventilation system. The clinical benefit is improved drug delivery.

STORAGE – HANDLING – DISPOSAL

- The recommended shelf life of MinimHal® is that recommended for the mechanical ventilation system with which it is associated, not to exceed 7 days, or less if visual examination indicates a need for replacement.
- MinimHal® should be discarded along with the mechanical ventilation system in accordance with institutional guidelines.
- MinimHal® should never be used on more than one circuit and should be disposed of in the same manner as the circuit on which it was installed.
- Any total or partial replacement of the circuit must lead to the replacement of MinimHal®.
- MinimHal® must be stored in a dry place.

WARNINGS – CONTRAINDICATIONS

- This device requires to be used by qualified staff to take care of patients who are under respiratory assistance.
- MinimHal® is a single patient product. Do not clean or sterilize. Do not reuse with other patients to avoid any risk of contamination.
- Any heat and moisture exchanger, located between the spacer and the patient, must be temporarily removed during the use of MinimHal®.
- Do not use MinimHal® with ventilation pressures exceeding 35 mbar, without permanent and careful attention.

The trademark Respimat® belongs to Boehringer Ingelheim. The trademark MinimHal® belongs to OptimHal-ProtectSom SAS. Patents belongs to Tours University and OptimHal-ProtectSom SAS together.

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