RC-FIT® OIMT

Mundstück / Mouthpiece / Boquilla / Boccaglio / Embout buccal

Inlay mit Ventil / Inlay with valve / Inserto con válvula / Inserto con valvola / Insert avec valve

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#### △ WARNING

Please clean before the first use! Please ensure that the RC-FIT® OIMT is being used in accordance with these instructions for use. Its incorrect use (off-label use) will result in an exclusion of liability for the manufacturer.

#### IMPORTANT INFORMATION

Before using the RC-FIT® OIMT for the first time, please read the Instructions for Use carefully. If you have additional questions regarding its usage, please consult with your qualified medical practitioner.

The RC-FIT® OIMT is well-suited for users starting at 3 years of age in home and clinical setting. The patient must be able to purposely coordinate inhalation and exhalation through the mouth. Children should only undergo training with the assistance of an adult. In the case of severely impaired patients, the training should be conducted together with a nurse or physiotherapist. Paraplegic patients or persons who suffer from neuromuscular diseases should only undergo training on the advice of qualified medical practitioners. Please abstain from any food intake shortly before and during the training.

Action to be taken in the event of a severe coughing fit: Following intake of your medication, try to work slowly with the respective RC-FIT® OIMT in its initial setting, in order to also dilate the bronchi mechanically.

Important: All individual parts must be stored out of the reach of children. Training following an eye operation (cataracts or glaucoma) should only start beginning the 4th week post-surgery.

After every use, please allow the respective RC-FIT® OIMT to dry as moisture may form in the product with long-term training.

#### USE OF OSCILLATING INHALATION TRAINING:

- several times a day and as needed
- approx. 10–30 breaths
- inhalation only through the mouth

The RC-FIT® OIMT is a sound-generating, oscillating, inspiratory trainer. This allows prevention, rehabilitation as well as the targeted improvement in fitness can be supported.

#### SCOPE OF DELIVERY (MATERIAL COMPOSITION)

##### RC-FIT® OIMT

Mouthpiece	MOUTHPIECE GREEN with 4 intensity levels [in mm: 6/5/4/3,75] (polypropylene)
Inlay with integrated inhalation valve	transparent (polypropylene/rust-free stainless steel)
Intensity regulator	grey (polypropylene)
Lid	grey (polypropylene)

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#### RC-FIT® OIMT

##### A) TRAINING WITH VIBRATIONS

The same principle applies to sports and everyday life: whoever has more oxygen available has a greater ability to function.  
The inlay with integrated inhalation valve produces oscillations that are acoustically perceived as sounds. This trains the strength and endurance of the respiratory muscles easily and efficiently, and loosens stubborn secretions.

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##### USE

Use of the RC-FIT® OIMT is not dependent on the force of gravity. The device may thus be used in a reclining, sitting or standing position. For the best effect, however, using the device while standing is preferable. Drink a glass of water before training in order to keep your airways moist.

- I. Please inspect the device for foreign bodies or damages prior to every use. If any parts are missing or damaged, replace them immediately.
- II. Always start your first training session in the lowest intensity level "MINIMUM" (largest possible triangle). The smaller the triangle becomes, the higher your training load (intensity). The RC-FIT® OIMT has a total of 4 intensity levels.

##### Average pressure during inhalation

intensity level	Flow rate (pressure difference to ambient pressure)			
10 l/min	15 l/min	20 l/min	25 l/min	

<tbl\_r cells="5" ix="5" maxcspan="

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The proprietor of the trademarks and registered trademarks is the company CEGLA Medizintechnik GmbH.