




RC-FIT® OIMT

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Please clean before the first use! Please ensure that the RC-FIT® OIMT is being used in accordance with these instructions for use. Its incorrect use (off-label use) will result in an exclusion of liability for the manufacturer.

IMPORTANT INFORMATION

Before using the RC-FIT® OIMT for the first time, please read the Instructions for Use carefully. If you have additional questions regarding its use, please consult with your qualified medical practitioner. The RC-FIT® OIMT is well-suited for users starting at 3 years of age in home and clinical setting. The patient must be able to purposely coordinate inhalation and exhalation through the mouth. Children should only undergo training with the assistance of an adult. In the case of severely impaired patients, the training should be conducted together with a nurse or physiotherapist. Paraplegic patients or persons who suffer from neuromuscular diseases should only undergo training on the advice of qualified medical practitioners. Please abstain from any food intake shortly before and during the training. Action to be taken in the event of a severe coughing fit: Following intake of your medication, try to work slowly with the respective RC-FIT® OIMT in its initial setting, in order to also dilate the bronchi mechanically. **Important: All individual parts must be stored out of the reach of children.** After every use, please allow the respective RC-FIT® OIMT to dry as moisture may form in the product with long-term training.

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- several times a day and as needed
- approx. 10–30 breaths
- inhalation only through the mouth

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The RC-FIT® OIMT is a sound-generating, oscillating, inspiratory trainer. This allows prevention, rehabilitation as well as the targeted improvement in fitness can be supported.

SCOPE OF DELIVERY (MATERIAL COMPOSITION)	
RC-FIT® OIMT	
Mouthpiece	MOUTHPIECE GREEN with 4 intensity levels [in mm: 6/5/4/3.75] (polypropylene)
Inlay with integrated inhalation valve	transparent (polypropylene/ rust-free stainless steel/Nickel silver)
Intensity regulator	Grey (polypropylene)
Lid	Grey (polypropylene)

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Training with the RC-FIT® OIMT produces a negative pressure in the lungs during inhalation. Therefore, contraindications for its use are:

- Congestive heart failure
- Untreated pneumothorax
- Respiratory tract infections
- Sinusitis
- In the case of a traumatic pneumothorax or rib fractures, these should be fully healed first

The following relative contraindications should be discussed with the doctor before use:

- untreated bronchial asthma
- untreated emphysema
- neuromuscular diseases

Note: Training following an eye operation (cataracts or glaucoma) should only start beginning the 4th week post-surgery.

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If inhalation training is too intense and results in the patient becoming winded, or experiencing vertigo or an excessively high pulse, please stop training. The next training unit should be restarted with a lower intensity level, or with a smaller number of breaths or shorter breath lengths.

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PURPOSE OF USE

The RC-FIT® OIMT trains the strength and endurance of the inhalation muscles, and loosens stubborn secretions in the case of acute and chronic diseases of the lower respiratory tract.

Furthermore, the device has several applications related to breathing and its psychophysiological impact, implying a healthier life and an increased well-being. At the core is the promotion of mindfulness and the stress reduction, the alleviation of pain or pain perception, support for mild psychological conditions such as depression and anxiety, and, as a result, a possible improvement in sleep.

Duration of use
If used daily, we recommend the exchange of the RC-FIT® after one year. Please notify the manufacturer of all serious incidents occurring in connection with the device.

Storage
Store the RC-FIT® in dry condition and protect the assistive device from sunlight and dust, among other things.

Disposal
The RC-FIT® can be disposed of with household waste, as long as not prescribed otherwise by the respective country's disposal regulations.

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RC-FIT® OIMT

A) TRAINING WITH VIBRATIONS

The same principle applies to sports and everyday life: whoever has more oxygen available has a greater ability to function. The device may thus be used in a reclining, sitting or standing position. For the best effect, however, using the device while standing is preferable. Drink a glass of water before training in order to keep your airways moist.

USE

Use of the RC-FIT® OIMT is not dependent on the force of gravity. The device may thus be used in a reclining, sitting or standing position. For the best effect, however, using the device while standing is preferable. Drink a glass of water before training in order to keep your airways moist.

I. Please inspect the device for foreign bodies or damages prior to every use. If any parts are missing or damaged, replace them immediately.

II. Always start your first training session in the lowest intensity level "MINIMUM" (largest-possible triangle). The smaller the triangle becomes, the higher your training load (intensity). The RC-FIT® OIMT has a total of 4 intensity levels.

Average pressure during inhalation

intensity level	Flow rate (pressure difference to ambient pressure)	10 l/min	15 l/min	20 l/min	25 l/min
I	-1 cmH ₂ O	-1,1 cmH ₂ O	-1,3 cmH ₂ O	-1,3 cmH ₂ O	-1,3 cmH ₂ O
II	-1,1 cmH ₂ O	-1,3 cmH ₂ O	-1,4 cmH ₂ O	-1,4 cmH ₂ O	-1,5 cmH ₂ O
III	-1,5 cmH ₂ O	-1,7 cmH ₂ O	-1,9 cmH ₂ O	-2,4 cmH ₂ O	-2,4 cmH ₂ O
IV	-3,2 cmH ₂ O	-4,3 cmH ₂ O	-5,7 cmH ₂ O	-10,5 cmH ₂ O	-10,5 cmH ₂ O

Frequency during inspiration: 165 Hz +/-5 %

III. A) Place the RC-FIT® OIMT in your mouth and close your lips around the mouthpiece. Please do not close the triangle with your fingers.

Breathe in through the mouth and out through the nose.

A) If the sound at the start of each breathing manoeuvre is loud and then quickly becomes quieter, adjust your breathing accordingly.

A) Make sure that the sound when breathing through the RC-FIT® OIMT is uniform and continuous. In such a case you will know that you have surmounted the training load.

B) Concentrate on the sounds occurring with your breathing. Try to produce these sounds to be as balanced as possible, for as long as possible.

A) Depending on your training condition, take a 1-minute break after 1½ minutes. One training unit should consist of 10-30 breaths, several times a day and as needed. In the case of Parkinson's with dysphagia, training should be conducted 5 x daily with 5 exhalations each time.

Is the training is too light, or the sound is too loud even after 15 breaths, increase the intensity (resistencia) by rotating the mouthpiece, see Fig. II.

CLEANING AND DISINFECTION (Boiling)

If used regularly, clean and disinfect the device once a week, please clean daily if there is a lot of mucus. If necessary (e.g. if no more sounds are produced), insert parts of the RC-FIT® OIMT in a warm, light detergent solution and then boil out.

1. Cleaning and disinfection

For cleaning purposes, pull off the mouthpiece from the intensity regulator. Then, carefully remove the inlay with integrated valve from the intensity regulator.

IMPORTANT: Try to avoid any contact with the integrated valve! Wash all parts in a solution of lukewarm clean water and standard liquid soap detergent. Rinse with clear water. After this cleaning, disinfect all parts in a vaporisator or in boiling water for approx. 5 minutes.

2. Drying

In order to prevent mould formation, all parts must be completely dried. To accelerate this process, shake off any excess water before setting the parts down. Allow all parts to air dry. Do not rub dry the inlay with towels.

3. Assembly

Before assembly, please wash your hands and ensure that all parts of the RC-FIT® OIMT are clean. There may be no dust residue or scattered particle deposits, as these could be breathed in and reach the lungs.

FREQUENCY OF USE

Recommended therapy intervals for the respiration: OIMT: several times a day and as needed, approx. 10-30 breaths each time

TIPS:

- Integrate the therapy and training into your daily routine, especially when you notice very little mucus, e.g. always before brushing your teeth.
- To obtain a better effect from your medication spray, integrate this inhalation into therapy with the RC-FIT® OIMT.

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The proprietor of the trademarks and registered trademarks is the company CEGLA Medizintechnik GmbH. We reserve the right to make technical and visual changes as well as to correct printing errors.