
RC-FIT® PEP/IMT

△ WARNHINWEISE

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△ WARNING

Please clean before the first use! Please ensure that the RC-FIT® PEP/IMT is being used in accordance with these Instructions for Use. Its incorrect use (off-label use) will result in an exclusion of liability for the manufacturer.

IMPORTANT INFORMATION

Before using the RC-FIT® PEP/IMT for the first time, please read the Instructions for Use carefully. If you have additional questions regarding its usage, please consult with your qualified medical practitioner. The RC-FIT® PEP/IMT is well-suited for users starting at 3 years of age in home and clinical setting. The patient must be able to purposely coordinate inhalation and exhalation through the mouth. Children should only undergo training/therapy with the assistance of an adult. In the case of severely impaired patients, the training/therapy should be conducted together with a nurse or physiotherapist. Paraplegic patients or persons who suffer from neuromuscular diseases should only undergo training/therapy on the advice of qualified medical practitioners. Please abstain from any food intake shortly before and during the training/therapy. Action to be taken in the event of a severe coughing fit: Following intake of your medication, try to work slowly with the respective RC-FIT® PEP/IMT in its initial setting, in order to also dilate the bronchi mechanically. **Important: All individual parts must be stored out of the reach of children.** After every use, please allow the respective RC-FIT® PEP/IMT to dry as moisture may form in the product with long-term training/therapy.

<b>USE OF PEP THERAPY:</b>	<b>USE OF IMT INHALATION TRAINING:</b>
• 2-3 x daily	• several times a day and as needed
• with prolonged exhalation for approx. 20 breaths each time	• approx. 1-3 minutes
• exhalation only through the mouth	• inhalation only through the mouth

For the PEP/IMT combination, please inhale and exhale through the mouth!

SCOPE OF DELIVERY (MATERIAL COMPOSITION)

RC-FIT® PEP/IMT	
Mouthpiece	MOUTHPIECE TRANSPARENT with 7 intensity levels [in mm: 4/3.75/3.5/3.25/3/2.5/2] (polypropylene)
Intensity regulator	Grey (polypropylene)
Lid	Grey (polypropylene)

△ CONTRAINDICATIONS

- Training and therapy with the RC-FIT® PEP/IMT produce in the lungs a negative pressure during inhalation and an overpressure during exhalation. Therefore, contraindications for its use are:
- Presence of pneumothorax, blood-expectorating cough (massive haemoptysis), and cranial pressure symptoms > 20 mmHg
- Decompensated cardiovascular diseases
- Acute thromboses
- Elevated coagulation factor
- Acute blood-expectorating cough (haemoptysis)
- Possible eardrum injuries
- Pulmonary barotrauma
- Significant cardiovascular insufficiency
- Massive restriction of the respiratory pump
- Skin emphysema
- Severe infection with fever

The following relative contraindications should be discussed with your treating healthcare professional:

- Tracheo-esophageal fistula
- Bronchial asthma
- Pulmonary emphysema
- Sinus infection
- Repeated, substantial nosebleeds
- Granuloma – growths in the mucous membranes near the tracheostoma
- after Surgery or injuries to the mouth, face, skull, or oesophagus

**Note:** Training/therapy following an eye operation (cataracts or glaucoma) should only start beginning the 4th week post-surgery.

SIDE EFFECTS

An exhalation pressure of more than 20 cm H<sub>2</sub>O may cause patients sensitive to increased transpulmonary pressure to experience one or more of the listed side effects:

- Hyperventilation (excessive increase in respiration) may occasionally cause slight vertigo or numbness in the mouth and hands. These symptoms are not dangerous and diminish once the respiration has normalised.
- Hypoventilation with hypercapnia (elevated CO<sub>2</sub> content).

**Other possible side effects**

- Elevated intracranial pressure
- Cardiovascular risk
- > Insufficient oxygen supply to the heart
- > Decreased venous return
- Swallowing of air with the possibility of vomiting and aspiration of vomit

If inhalation training is too intense and results in the patient becoming winded, or experiencing vertigo or an excessively high pulse, please stop training. The next training unit should be restarted with a lower intensity level, or with a smaller number of breaths or shorter breath lengths.

PURPOSE OF USE

The RC-FIT® PEP/IMT trains the strength and endurance of the inhalation and exhalation muscles, and keeps the airways open in the case of acute and chronic diseases of the lower respiratory tract. The RC-FIT® PEP/IMT can effectively train the pursed lips breathing technique. This training of the respiratory tract holds the airways open during exhalation using resistance. In this manner, a bronchial collapse is prevented. Additionally, the RC-FIT® PEP/IMT is a training device and an aid in emergency cases. Furthermore, the device has several applications related to breathing and its psychophysiological impact, implying a healthier life and an increased well-being. At the core is the promotion of mindfulness and stress reduction, the alleviation of pain or pain perception, support for mild psychological conditions such as depression and anxiety, and, as a result, a possible improvement in sleep.

**Duration of use**  
If used daily, we recommend the exchange of the RC-FIT® PEP/IMT after one year. Please notify the manufacturer of all serious incidents occurring in connection with the device.  
**Storage**  
Store the RC-FIT® PEP/IMT in dry condition and protect the assistive device from sunlight and dust, among other things.  
**Disposal**  
The RC-FIT® PEP/IMT can be disposed of with household waste, as long as not prescribed otherwise by the respective country's disposal regulations.

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RC-FIT® PEP/IMT

A) TRAINING & THERAPY

The same principle applies to sports and everyday life: whoever has more oxygen available has a greater ability to function.  
By only breathing in through the RC-FIT® PEP/IMT, the respiratory muscles can be trained. By also breathing out through the RC-FIT® PEP/IMT, the airways can be expanded.

B) ANCHOR FOR SHORTNESS OF BREATH

As with the pursed lips technique, the RC-FIT® PEP/IMT can stabilise the bronchi in cases of shortness of breath, through a deliberately slow exhalation. Hyperventilation can be reduced or prevented. Thanks to its practical trouser pocket shape, the device is always close at hand in the lowest intensity level.  
Practice removing with the RC-FIT® PEP/IMT and the Use pursed lip technique about 60 times to be optimal for emergencies to be armed.

USE

Use of the RC-FIT® PEP/IMT is not dependent on the force of gravity. The device may thus be used in a reclining, sitting or standing position. For the best effect, however, using the device while standing is preferable. Drink a glass of water before training in order to keep your airways moist.

I. Please inspect the device for foreign bodies or damages prior to every use. If any parts are missing or damaged, replace them immediately.

II. Always start your first training session in the lowest intensity level “MINIMUM” (largest-possible triangle).  
The smaller the triangle becomes, the higher your training and therapy load (intensity). The RC-FIT® PEP/IMT has a total of 7 intensity levels.

Average pressure during exhalation/ inhalation

intensity level	Flow rate (pressure difference to ambient pressure)			
	10/-10 l/min	15/-15 l/min	20/-20 l/min	25/-25 l/min
I	0,9/-1,1 cmH <sub>2</sub> O	2/-2,3 cmH <sub>2</sub> O	3,6/-4,1 cmH <sub>2</sub> O	5,5/-6,4 cmH <sub>2</sub> O
II	1,1/-1,3 cmH <sub>2</sub> O	2,5/-2,7 cmH <sub>2</sub> O	4,4/-4,8 cmH <sub>2</sub> O	6,6/-7,5 cmH <sub>2</sub> O
III	1,4/-1,5 cmH <sub>2</sub> O	2,9/-3,4 cmH <sub>2</sub> O	5/-5,8 cmH <sub>2</sub> O	7,7/-8,9 cmH <sub>2</sub> O
IV	1,9/-1,9 cmH <sub>2</sub> O	4,1/-4,2 cmH <sub>2</sub> O	7,1/-7,2 cmH <sub>2</sub> O	10,9/-11,1 cmH <sub>2</sub> O
V	2,2/-2,2 cmH <sub>2</sub> O	4,6/-4,6 cmH <sub>2</sub> O	8,1/-8,1 cmH <sub>2</sub> O	12,3/-12,3 cmH <sub>2</sub> O
VI	3,2/-3,2 cmH <sub>2</sub> O	7,2/-6,9 cmH <sub>2</sub> O	12,3/-11,7 cmH <sub>2</sub> O	18,7/-17,9 cmH <sub>2</sub> O
VII	5,3/-5,6 cmH <sub>2</sub> O	11,5/-11,8 cmH <sub>2</sub> O	19,4/-20,3 cmH <sub>2</sub> O	27,7/-30,1 cmH <sub>2</sub> O

III. A) Place the RC-FIT® PEP/IMT in your mouth and close your lips around the mouthpiece. Please do not close the triangle with your fingers.

PEP: Inhalation through the nose and exhalation only through the mouth  
IMT: Inhalation only through the mouth and exhalation through the nose or inhalation and exhalation via the RC-FIT®

A) Depending on your training condition, take a 1-minute break after 1½ minutes. One training unit should consist of 10–30 breaths 3 x daily. In the case of Parkinson's with dysphagia, training of 5 x daily with 5 exhalations each time should be conducted.

C) To train for cases of shortness of breath or while experiencing shortness of breath always use only the “MINIMUM” intensity level. Try to breathe in the goalkeeper position or the cart driver position calmly and in a prolonged manner, via the RC-FIT® PEP/IMT. Repeat this breathing cycle approx. 10 times if possible.

If training feels too easy, increase the intensity (resistance) by rotating the mouthpiece, see Fig. II.

CLEANING AND DISINFECTION (Boiling)

If used regularly, clean and disinfect the device once a week, please clean daily if there is a lot of mucus.

1. Cleaning and disinfection

For cleaning purposes, pull off the mouthpiece from the intensity regulator.  
Wash all parts in a solution of lukewarm clean water and standard liquid soap detergent. Rinse with clear water. After this cleaning, disinfect all parts in a vaporisator or in boiling water for approx. 5 minutes.

2. Drying

In order to prevent mould formation, all parts must be completely dried. To accelerate this process, shake off any excess water before setting the parts down. Allow all parts to air dry.

3. Assembly

Before assembly, please wash your hands and ensure that all parts of the RC-FIT® PEP/IMT are clean. There may be no dust residue or scattered particle deposits, as these could be breathed in and reach the lungs.

FREQUENCY OF USE

Recommended therapy intervals for the respiration:  
PEP: 2–3 x daily with prolonged exhalation for approx. 20 breaths each time  
IMT: several times a day and as needed, approx. 1–3 minutes

TIPS:

- ▶ Integrate the therapy and training into your daily routine, especially when you notice very little mucus, e.g. always before brushing your teeth.
- ▶ To obtain a better effect from your medication spray, integrate this inhalation into therapy with the RC-FIT® PEP/IMT.

ACCESSORIES

The RC-Mask Tracheostomy allows tracheostomised patients to undergo respiratory therapy and respiratory training with the RC-FIT® PEP/IMT. The RC-Mask Tracheostomy is simply placed over the mouthpiece of the RC-FIT® model and then connected with the open tracheostomy, i.e. the cannula.

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