Chambre d'inhalation pour circuit de ventilation mécanique.
Inhalation chamber for mechanical ventilation circuit.
Inhalationskammer für den maschinellen Beatmungskreislauf.
Camera di inalazione per circuito di ventilazione meccanica.

PI-NT-COM-0003 / 20230105

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INTENDED USE

CombiHaler® is an inhalation chamber intended for adult patient with invasive or non-invasive mechanical ventilation who have been prescribed aerosolized treatment from a pressurized metered-dose inhaler and/or a vibrating mesh nebulizer.

NOTES

- CombiHaler® should be used by clinicians qualified in the care of mechanically ventilated patients.
- Please read the instructions for use of CombiHaler®, the metered dose inhaler and the nebulizer and keep them available at all times of use.

COMPOSITION (FIG.1)

a. Male connector 22 mm
b. Female connector 22 mm
c. One opening for MDI
d. One opening for vibrating mesh nebulizer
e. Two sealed caps
f. One nozzle
g. A RespiMat® inhaler adapter can also be supplied.

REQUIREMENTS BEFORE USE

When putting the device into operation, it is imperative to ensure that:

- The device does not contain any foreign body, to avoid any risk of damage to the patient's airways or obstruction of the medical device.
- The caps are correctly positioned when MDI and/or nebulizer is not used.

INSTRUCTIONS FOR USE

Checking the operation of the nozzle
Remove the MDI canister from its actuator. Open the orifice cap (c) and insert the canister (Fig.2).
Operate the MDI and check that the nozzle correctly ensures its function (visually test with a puff of drugs) (Fig.3).

Insertion in the circuit
In a traditional adult adult dual limb of mechanical ventilation circuit, the CombiHaler® must be inserted on the inspiratory limb. The arrow on the device must point towards the patient (FIG. 4). If, for technical reason, CombiHaler® would be inserted between the Y-piece and the patient, it's imperative that the inhalation chamber is removed from the system immediately after taking the medication and is never left in place while the circuit is being used.

DRUG ADMINISTRATION BY PRESSURIZED METERED DOSE INHALER

Step 1:
Shake the MDI canister before use in accordance with the instructions provided by the manufacturer. Remove the cap from the opening (c) and insert the canister (Fig.5).

Step 2:
Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the MDI (Fig.6).

Step 3:
Remove the canister from the MDI port and replace the cap on the orifice, shake the canister and repeat steps 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the aerosol dispenser for the waiting time between each dose. At the end of the drug administration, remove the canister from the opening MDI and replace the cap on the orifice (c) (Fig.7).

DRUG ADMINISTRATION BY VIBRATING MESH NEBULIZER

Remove the cap from the orifice (d) and insert the vibrating mesh nebulizer (Fig.8). To avoid an opening of the circuit, this step could be done before the insertion of CombiHaler® into the circuit. Pour the drug solution into the nebulizer reservoir. For proper operation of the nebulizer, follow the instructions provided by the manufacturer.

Aerogen® vibrating mesh nebulizers may remain on the circuit as recommended by the manufacturer.

DRUG ADMINISTRATION BY VIBRATING MESH NEBULIZER

Remove the cap from the orifice (d) and insert the vibrating mesh nebulizer (Fig.8). To avoid an opening of the circuit, this step could be done before the insertion of CombiHaler® into the circuit. Pour the drug solution into the nebulizer reservoir. For proper operation of the nebulizer, follow the instructions provided by the manufacturer.

Aerogen® vibrating mesh nebulizers may remain on the circuit as recommended by the manufacturer.

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DRUG ADMINISTRATION WITH THE RESPIMAT® INHALER

Step 1:

Remove the cap from the orifice (d) and insert the adapter Respimat®. Insert the Respimat® inhaler into the adapter (Fig.9).

Step 2:

To operate the inhaler, follow the instructions provided by the manufacturer. Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the inhaler.

Step 3:

Repeat step 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the inhaler dispenser for the waiting time between each dose. At the end of the drug administration, remove the inhaler and the adapter from the opening (d) and replace the cap on the orifice (d).

CLINICAL BENEFITS – DEVICE PERFORMANCE

The chamber volume allows the propellant gas of the larger particles size to evaporate, producing smaller aerosol particles of the drug. Smaller particles improve the efficiency of aerosol delivery to the patient through the ventilation system. The clinical benefit is improved drug delivery.

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STORAGE – HANDLING – DISPOSAL

- The recommended shelf life of CombiHaler® is that recommended for the mechanical ventilation system with which it is associated, not to exceed 7 days, or less if visual examination indicates a need for replacement.
- CombiHaler® should be discarded along with the mechanical ventilation system in accordance with institutional guidelines.
- CombiHaler® should never be used on more than one circuit and should be disposed of in the same manner as the circuit on which it was installed.
- Any total or partial replacement of the circuit must lead to the replacement of CombiHaler®.
- CombiHaler® must be stored in a dry place.

WARNINGS – CONTRAINDICATIONS

- This device requires to be used by qualified staff to take care of patients who are under respiratory assistance.
- CombiHaler® is a single patient product. Do not clean or sterilize. Do not reuse with other patients to avoid any risk of contamination.
- To ensure the holding of the device, CombiHaler® must not be dislocated. In case of dislocation, CombiHaler® must not be used, and must be replaced.
- CombiHaler® is design to be inserted into an adult mechanical ventilation circuit.
- Any heat and moisture exchanger, located between the spacer and the patient, must be temporarily removed during the use of CombiHaler®.
- Do not use CombiHaler® with ventilation pressures exceeding 35 mbar, without permanent and careful attention.

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CombiHaler® is a trademark of OptimHal-ProtecSom SAS. The French and international patents that apply to CombiHaler® are the joint property of OptimHal-ProtecSom SAS and the university of Tours. Aerogen® is a trademark of Aerogen Compagny Limited. Respimat® is a trademark of Boehringer Ingelheim.

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