RC-FIT® CLASSIC

Mundstück / Mouthpiece / Boquilla /
Boccaglio / Embout buccalInlay mit integrierten
Ventilen / Inlay with valve /
Inserto con válvula / Inserto con valvola /
Insert avec valveIntensitätsregler / Intensity regulator /
Regulador de intensidad / Regolatore
di intensità / Régulateur d'intensité

Deckel / Lid / Tapa / Coperchio / Cache

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▲ WARNING

Please clean before the first use! Please ensure that the RC-FIT® CLASSIC is being used in accordance with these Instructions for Use. Its incorrect use (off-label use) will result in an exclusion of liability for the manufacturer.

IMPORTANT INFORMATION

Before using the RC-FIT® CLASSIC for the first time, please read the Instructions for Use carefully. If you have additional questions regarding its usage, please consult with your qualified medical practitioner. The RC-FIT® CLASSIC is well-suited for users starting at 3 years of age in home and clinical setting. The patient must be able to purposely coordinate inhalation and exhalation through the mouth. Children should only undergo training/therapy with the assistance of an adult. In the case of severely impaired patients, the training/therapy should be conducted together with a nurse or physiotherapist. Paraplegic patients or persons who suffer from neuromuscular diseases should only undergo training/therapy on the advice of qualified medical practitioners. Please abstain from any food intake shortly before and during the training/therapy. Action to be taken in the event of a severe coughing fit: Following intake of your medication, try to work slowly with the respective RC-FIT® CLASSIC in its initial setting, in order to also dilate the bronchi mechanically. **Important: All individual parts must be stored out of the reach of children.** Training/therapy following an eye operation (cataracts or glaucoma) should only start beginning the 4th week post-surgery. After every use, please allow the respective RC-FIT® Classic to dry as moisture may form in the product with long-term training/therapy.

USE OF OSCILLATING EXHALATION THERAPY AND INHALATION TRAINING

- At least 3 times a day (more if necessary)
- with prolonged exhalation for approx. 10–30 breaths each time
- inhalation and exhalation only through the mouth

SCOPE OF DELIVERY (MATERIAL COMPOSITION)

RC-FIT® CLASSIC	
Mouthpiece	MOUTHPIECE BLUE with 4 intensity levels [in mm: 6/5/4/3,75] (polypropylene)
Inlay with integrated valves	Transparent (polypropylene/rust-free stainless steel)
Intensity regulator	Grey (polypropylene)
Lid	Grey (polypropylene)

▲ CONTRAINDICATIONS

Training and therapy with the RC-FIT® CLASSIC produce in the lungs a negative pressure during inhalation, and an overpressure during exhalation. Therefore, contraindications for its use are:

- Respiratory diseases with fever
- Presence of pneumothorax, blood-expectorating cough (massive haemoptysis), and cranial pressure symptoms > 20 mmHg.
- Decompensated cardiovascular diseases
- Acute thromboses/elevated coagulation status during or after Covid-19 disease
- Acute blood-expectorating cough (haemoptysis)
- Surgery or injuries to the mouth, face, skull, or oesophagus
- Possible eardrum injuries
- Pulmonary barotrauma
- Significant cardiovascular insufficiency
- Massive restriction of the respiratory pump
- Repeated, substantial nosebleeds
- Skin emphysema
- Granuloma - growths in the mucous membranes near the tracheostoma

The following relative contraindications should be discussed with your treating healthcare professional:

- Tracheoesophageal fistula
- Bronchial asthma
- Pulmonary emphysema
- Sinus infection

SIDE EFFECTS

An exhalation pressure of more than 20 cm H₂O can cause patients sensitive to increased transpulmonary pressure to experience one or more of the listed side effects:

- Hyperventilation (excessive increase in respiration) may occasionally cause slight vertigo or numbness in the mouth and hands. These symptoms are not dangerous and diminish once the respiration has normalised.
- Hypoventilation with hypercapnia (elevated CO₂ content).

Other possible side effects

- Elevated intracranial pressure
- Cardiovascular risk
 - > Insufficient oxygen supply to the heart
 - > Decreased venous return
- Swallowing of air with increased probability of vomiting and aspiration of vomit

If inhalation training is too intense and results in the patient becoming winded, or experiencing vertigo or an excessively high pulse, please stop training. The next training unit should be restarted with a lower intensity level, or with a smaller number of breaths or shorter breath lengths.

PURPOSE OF USE

The RC-FIT® CLASSIC is intended for secretion mobilisation in the case of acute and chronic diseases of the lower respiratory tract. The RC-FIT® CLASSIC is a respiratory trainer and respiratory therapy device in one and help in an emergency.

Duration of use

If used daily, we recommend the exchange of the RC-FIT® CLASSIC after one year. Please notify the manufacturer of all serious incidents occurring in connection with the device.

Storage

Only store the RC-FIT® CLASSIC in dry condition and with closed lid. In this manner, the assistive device is protected from sunlight and dust, among other things.

Disposal

The RC-FIT® CLASSIC can be disposed of with household waste, as long as not prescribed otherwise by the respective country's disposal regulations.

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RC-FIT® CLASSIC

A) TRAINING & THERAPY WITH VIBRATIONS

The same principle applies to sports and everyday life: whoever has more oxygen available has a greater ability to function. The inlay with integrated valve produces oscillations that are acoustically perceived as sounds. These expand in the airways. As a result, shortness of breath is decreased during physical exertion. Mucilage and coughing in acute and chronically ill patients can be reduced. Through stabilisation of the airways more oxygen can be taken in. In addition, body function perception is achieved. The goal is to support the cleaning function of the airways in the best possible manner. This can be very effective within the scope of simple colds all they way to chronic pulmonary diseases. Aside from moistening the airways with fresh, moist air (if there is no fever present), inhalations, or moistening nasal sprays in order to keep the cleaning mechanism working, the PEP therapy device can also free the airways of stubborn mucus. The breeding ground for infections is thus minimised.

B) ANCHOR FOR SHORTNESS OF BREATH

As with the pursed lips technique, the RC-FIT® CLASSIC can stabilise the bronchi in cases of shortness of breath, through a deliberately slow exhalation. Hyperventilation can be reduced or prevented. Thanks to its practical trouser pocket shape, the device is always close at hand in the lowest intensity level. Practice removing with the RC-FIT® CLASSIC and the Use pursed lip technique about 60 times to be optimal for emergencies to be armed.

USE

Use of the RC-FIT® CLASSIC is not dependent on the force of gravity. The device may thus be used in a reclining, sitting, standing or moving position. For the best effect, however, using the device while standing is preferable. Drink a glass of water before training in order to keep your airways moist.

- Please inspect the device for foreign bodies or damages prior to every use. If any parts are missing or damaged, replace them immediately.
- Always start your first training session in the lowest intensity level "MINIMUM" (largest-possible triangle). The smaller the triangle becomes, the higher your training and therapy load (intensity). The RC-FIT® CLASSIC has a total of 4 intensity levels.

intensity level	Average pressure during exhalation/ inhalation			
	Flow rate (pressure difference to ambient pressure)			
	10/-10 l/min	15/-15 l/min	20/-20 l/min	25/-25 l/min
I	1,7/-1,4 cmH ₂ O	3,3/-3,2 cmH ₂ O	5,4/-4,2 cmH ₂ O	8,4/-5,1 cmH ₂ O
II	1,9/-1,7 cmH ₂ O	3,9/-3,7 cmH ₂ O	6,6/-5,3 cmH ₂ O	10/-6 cmH ₂ O
III	2,2/-2,2 cmH ₂ O	4,9/-4,8 cmH ₂ O	8,4/-8,3 cmH ₂ O	12,6/-12,1 cmH ₂ O
IV	2,2/-2,3 cmH ₂ O	5/-5,3 cmH ₂ O	8,4/-8,7 cmH ₂ O	12,9/-13,8 cmH ₂ O

Frequency during inspiration: 165 Hz+/5%
Frequency during exhalation: 130 Hz+/5%

III. A) Place the RC-FIT® CLASSIC in your mouth and close your lips around the mouthpiece. Please do not close the triangle with your fingers.

Breathe in and out through the mouth.

A) If the sound at the start of each breathing manoeuvre is loud and then quickly becomes quieter, please adapt your breathing accordingly.

A) Make sure that the sound when breathing through the RC-FIT® CLASSIC is uniform and continuous. In such a case you will know that you have surmounted the training load. Concentrate on the sounds occurring with your breathing. Try to produce these sounds to be as balanced as possible, for as long as possible.

A) Depending on your training condition, take a 1-minute break after 1½ minutes. One training unit should consist of 10–30 breaths, 3 x daily. In the case of Parkinson's with dysphagia, training of 5 x daily with 5 exhalations each time should be conducted.

B) To train for cases of shortness of breath or while experiencing shortness of breath always use only the "MINIMUM" intensity level. Try to breathe in the goalkeeper position or the car driver position calmly and in a prolonged manner, via the RC-FIT® CLASSIC. Repeat this breathing cycle approx. 10 times if possible.

If the training is too light, or the sound is too loud even after 30 breaths, increase the intensity (resistance) by rotating the mouthpiece, see Fig. II.

CLEANING

1. **Cleaning**
For cleaning purposes, pull off the mouthpiece from the intensity regulator. Then, carefully remove the inlay with integrated valve from the intensity regulator. IMPORTANT: Try to avoid any contact with the integrated valve! Wash all parts in a solution of lukewarm clean water and standard liquid soap detergent. Rinse with clear water. After this cleaning, all parts may be cleaned in a vaporiser or for approx. 5 minutes in boiling water. If used regularly, clean the device twice a week. Thoroughly once a week, please clean daily if there is a lot of mucus. If necessary (e.g. if no more sounds are produced), insert parts of the RC-FIT® CLASSIC in a warm, light detergent solution and then boil out.

2. **Drying**
In order to prevent mould formation, all parts must be completely dried. To accelerate this process, shake off any excess water before setting the parts down. Allow all parts to air dry. Do not rub dry the inlay with towels.

3. **Assembly**
Before assembly, please wash your hands and ensure that all parts of the RC-FIT® CLASSIC are clean. There may be no dust residue or scattered particle deposits, as these could be breathed in and reach the lungs.

FREQUENCY OF USE

Recommended therapy intervals for the respiration:
RC-FIT® CLASSIC: 3 x daily 10–30 breaths, 7 days a week. Additional usage when you notice heavy mucus obstruction or occasional shortness of breath.

TIPS:
▶ Integrate the therapy and training into your daily routine, especially when you notice very little mucus, e.g. always before brushing your teeth.
▶ To obtain a better effect from your medication spray, integrate this inhalation into therapy with the RC-FIT® CLASSIC.

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The proprietor of the trademarks and registered trademarks is the company CEGLA Medizintechnik GmbH. We reserve the right to make technical and visual changes as well as to correct printing errors.