MinimHal®


EN Spacer with low volume (60 mL) for mechanical ventilation circuit

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COMPOSITION

- Cone connector 15 mm
- Socket connector 15 mm
- One opening for MDI
- One opening for vibrating mesh nebulizer
- Two sealed caps
- Nozzle
- A RespiMat® inhaler adapter can also be supplied

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Description - Indications

MinimHal® is an inhalation chamber intended for invasive and non-invasive mechanical ventilation circuits. This device was designed for the administration of aerosol drugs for mechanically ventilated patients. It allows the use of both Metered Dose Inhalers (MDI) and vibrating mesh nebulizers. MinimHal® does not weight down on the ventilation circuit and allows monitoring of medication administration.

Indications for use

- The bag must be opened only when you intend to use MinimHal®. Once out of its packaging, this product can be used immediately.
- When starting, it is necessary to ensure that :
 - The device does not contain any foreign objects.
 - The MDI nozzle correctly ensures its function (visually test with a puff of drugs).
- During use, it is necessary to ensure that the caps are correctly positioned when MDI and/or nebulizer is not used.

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



Instructions for use

Checking the operation of the nozzle

Remove the MDI canister from its actuator. Open the orifice cap (c) and insert the canister.

Operate the MDI and inspect the nozzle for proper function. Look at the medication.

Insertion in the circuit

In a traditional adult (with two adapters) or pediatric dual-limb mechanical ventilation circuit, insert MinimHal® on the inspiratory branch of the circuit or after the Y piece*. The cone connector (a) must be plugged to the inspiratory branch (ill.1) or to the Y piece (ill.2). The socket connector (b) must be plugged to the Y piece (ill.1) or to the patient interface.

**When using MinimHal® after the Y piece, make sure that the dead volume (to which the chamber's volume is added - 60 mL) is compatible with the tidal volume delivered. At the end of the administration of medication, disconnect MinimHal® from the circuit.*

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Drug Administration by pressurized dose metered





Shake the MDI canister before use in accordance with the instructions provided by the manufacturer. Remove the cap from the opening (c) and insert the canister.

Hold the chamber in your hand in a perfectly horizontal position. At the beginning of the inspiration, operate the MDI.

Remove the canister from the MDI port and replace the cap on the orifice, shake the canister and repeat steps 1 and 2 according to the number of doses required. Follow the instructions of the manufacturer of the aerosol dispenser for the waiting time between each dose. At the end of the drug administration, remove the canister from the opening MDI and replace the cap on the orifice (c).

Administration of drugs by vibrating mesh nebulizer

Remove the cap from the orifice (d) and insert the vibrating mesh nebulizer. Pour the drug solution into the nebulizer reservoir. For proper operation of the nebulizer, follow the instructions provided by the manufacturer.

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


As soon as the medication is complete, remove the nebulizer and put the cap back on the orifice (d).

Drug Administration with the Respimat® inhaler

Remove the cap from the orifice (d) and insert the adapter Respimat®. Insert the Respimat® inhaler into the adapter. To operate the inhaler, follow the instructions provided by the manufacturer.

Repeat step 1 according to the number of doses required. Follow the instructions of the manufacturer of the MDI for the waiting time between each dose. As soon as the medication is complete, remove the inhaler and adapter and put the cap back on the orifice (d).

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Aerodynamics particle size characteristics (Results are expressed as means ± SD)

	With vibrating mesh nebulizer	With pMDI
Mass median aerodynamic diameter (MMAD) (µm)	1.86 ± 0.55	2.05 ± 0.11
Fines particles dose (< 5 µm) (µg)	2174 ± 418	252 ± 46
Geometric standard deviation	4.15 ± 1.24	1.99 ± 0.12

The MinimHal® spacer was evaluated for particle size characteristics using a new generation cascade impactor, at a constant flow rate of 28.3L/min. MinimHal® was connected to the cascade impactor through a 7.5 mm inside diameter endotracheal tube. MinimHal®'s performances were evaluated with a pMDI (Ventolin®, GSK; 10 doses of 100 µg were actuated in the device) and with the vibrating mesh nebulizer Aeroneb Solo® (Salbutamol Mylan, 2.5mg/2.5mL, nominal dose : 5000 µg).

Storage – handling – disposal

- MinimHal® is a single patient use product. Do not clean and reuse or sterilize.
- The recommended lifespan of MinimHal® is equivalent to the recommended service life for the mechanical ventilation circuits to which it is associated, without exceeding 7 days, or sooner if a visual examination indicates that it is necessary to make its change.
- MinimHal® must be disposed of at the same time as the mechanical ventilation circuit, in accordance with the instructions of the facility.
- MinimHal® must not be used on more than one circuit and must be destroyed under the same conditions as the circuit on which it was mounted.
- Any total or partial replacement of the circuit must result in the replacement of MinimHal®.
- Store away from moisture.

Warnings – contraindications

- This device requires to be used by qualified staff to take care of patients who are under respiratory assistance.
- Do not use with a pressure over 35 mbar without careful and permanent attention.
- Any heat and moisture exchanger, located between the spacer and the patient, must be temporarily removed during the use of MinimHal®.

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The trademark Respimat® belongs to Boehringer Ingelheim.
The trademark MinimHal® belongs to OptimHal-ProtectSom SAS.
OptimHal-ProtectSom SAS
MinimHal® est une marque de commerce d'OptimHal-ProtectSom SAS.
MinimHal® ist ein Warenzeichen von OptimHal-ProtectSom SAS.
MinimHal® è un marchio commerciale di OptimHal-ProtectSom SAS.